HI-FLOW™ SUCTION CANISTERS


(800, 1200, 2000 & 3000CC)


HYDROPHOBIC SUCTION CANISTERS

(800, 1200 & 2000CC)


SUCTION CANISTER	4	RECIPIENTE PER ASPIRAZIONE	16
АСПИРАЦИОНЕН КОНТЕЙНЕР	5	SUCTION CANISTER	17
SPREMNIK ZA USISAVANJE	6	IESŪCES TVERTNE	18
ODSÁVACÍ NÁDOBA	7	ĮSIURBIMO INDAS	19
SUGEBEHOLDER	8	KANISTER TA' ĠBID	20
OPVANGRESERVOIR VOOR VOCHT	9	POJEMNIK SSĄCY	21
ASPIRAATORI KOGUMISANUM	10	RESERVATÓRIO DE SUÇÃO	22
IMUSÄILIÖ	11	PUNGĂ PENTRU ASPIRARE	23
CARTOUCHE D'ASPIRATION	12	SACIA NÁDOBA	24
SAUGKANISTER	13	SESALNA POSODA	25
KANISTPO ANAPPOΦΗΣΗΣ	14	DEPÓSITO DE ASPIRACIÓN	26
SZÍVÓTARTÁLY	15	SUGKOPP	27


**SMALL
HYDROPHOBIC COVER
(800CC)**

**SMALL
HI-FLOW COVER
(800CC, 1200CC SOLO)**

**LARGE
HI-FLOW COVER
(1200CC, 2000CC, 3000CC)**

**LARGE
HYDROPHOBIC COVER
(1200CC, 2000CC)**

	EN	BG	HR	CZ	DA	NL
1	PATIENT	ПАЦИЕНТ	PACIJENT	PACIENT	PATIENT	PATIËNT
2	VACUUM	ВАКУУМ	VAKUUM	PODTLAK	VAKUUM	VACUÛM
3	TANDEM	ДВОЕН	TANDEM	TANDEM	TANDEM	TANDEM
4	POUR SPOUT	ДОЗАТОР	NASTAVAK ZA IZLIJEVANJE	VYLÉVACÍ OTVOR	TAPSNUDE	SCHENKTUIT

	ET	FI	FR	DE	EL	HU
1	PATSIENT	POTILAS	PATIENT	PATIENT	ΑΣΘΕΝΗΣ	BETEG
2	VAAKUM	VAKUUMI	VIDE	VAKUUM	KENO	VÁKUUM
3	TANDEMAVA	KAKSOIS	TANDEM	TANDEM	ΠΑΡΑΛΛΗΛΗ ΘΥΡΑ	TANDEM
4	VALAMISTILA	KAATONOKKA	BEC VERSEUR	AUSGIESSER	ΣΤΟΜΙΟ ΡΟΗΣ	KIÖNTÖNYÍLÁS

	IT	GA	LV	LT	MT	PL
1	PAZIENTE	PATIENT	PACIENTS	PACIENTAS	PAZJENT	PACJENT
2	VUOTO	VACUUM	VAKUUMS	VAKUUMAS	VAKWU	PODCIŚNIENIE
3	TANDEM	TANDEM	TANDĒMS	TANDEMAS	TANDEM	TANDEM
4	BECCUCCIO DOSATORE	POUR SPOUT	LEJAMTEKNE	IŠPYLIMO ANGA	ŽENNUNA GHATTFERRIGH	WLEW

	PT	RO	SK	SL	ES	SV
1	DOENTE	PACIENT	PACIENT	PACIENT	PACIENTE	PATIENT
2	VÁCUO	VACUUM	VÁKUUM	VAKUUM	VACÍO	VAKUUM
3	TANDEM	TANDEM	TANDEM	TANDEM	TÁNDEM	TANDEM
4	BOCAL	PORT DE GOLIRE	DÁVKOVACIA HUBICA	DULEC ZA IZLIVANJE	BOQUILLA DE VERTIDO	HÄLLTAPP

SUCTION CANISTER

INTENDED USE: *Intended for use with suction equipment to allow for the collection and disposal of bodily fluids or other materials removed during surgical or medical procedures.*

DIRECTIONS FOR USE

1. **Hi-Flow models only:** Gently shake lid to verify shut-off mechanism moves freely.
2. Place lid on canister and press firmly around entire perimeter.
3. Take right angle connectors (elbows) out of bags found in lid packages and apply tightly to patient and vacuum ports (except for 800cc and 1200cc Solo models).
Note: On Hi-Flow models only, the patient tube can be placed directly on patient port without an elbow if desired.
4. Apply pour spout cap firmly over pour spout and

apply tandem port cap tightly over tandem port if not being used (no tandem ports on 800cc and 1200cc Solo models).

5. Apply patient tubing to patient port and vacuum tube to vacuum port. Be sure tubing fits snugly.
6. Check all caps and connections for proper seal. Test the assembly for vacuum leaks by turning on vacuum source and occluding the patient tubing with finger or thumb.
7. To prevent fluid contamination in vacuum line, make sure vacuum line is attached to vacuum port.

DISPOSAL

1. Turn off vacuum source and disconnect all tubing.
2. Seal all ports with attached port caps.
3. Remove canister from bracket and transport to disposal area. **Do not lift canister by lid. The weight of the contents may cause the lid to separate from the canister.**
4. Dispose of according to your facility's standard operating procedures for biohazard disposal and both federal and local disposal regulations.

TROUBLESHOOTING

- Loss of Vacuum: Check that vacuum is on, canister is properly sealed, and that all connections are tight and tubing is not kinked. If loss continues, replace canister.

CAUTION

- Intended User: To be used by professionally trained healthcare personnel.
-  Single Use Only. Do not attempt to clean, sterilize or reuse canister. Possible consequences of reuse include: 1) implosion, 2) fluid bypass, and 3) exposure to pathogens.
- Do not exceed vacuum level of
 - 23" Hg (77.9 kPa) for 800cc Hydrophobic
 - 25" Hg (84.7 kPa) for all Hi-Flow products and 1200/2000cc Hydrophobic
-  Do not use for liposuction procedures.
- Do not apply continuous vacuum longer than 24 hours.
- Bemis components not designed to be compatible with other manufacturer's canisters or lids.
-  Not intended as a measuring device, only for general reference, not specific measurement.
-  Canister contents are considered potentially hazardous. Use appropriate PPE and handle accordingly.
-  Check expiration date on canister for use-by-date (if applicable).
-  Store in a dark place. Long term exposure to light may compromise product performance and result in breakage during use.

Any serious incident related to the use of this product, should be reported to both the manufacturer and the health authority/competent authority where the product is established.

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UK CA Responsible Person: Bemis Limited